
**Motiva®** Channel Separator  
Instructions for Use

**MOTIVA® CHANNEL SEPARATOR  
INSTRUCTIONS**

|                     |                               |
|---------------------|-------------------------------|
| <b>Ref:</b>         | MCS-20                        |
| <b>Description:</b> | Shaft 20 mm wide Length 310mm |
| <b>UDI-DI</b>       | 07350093650624                |

Figure 1. Motiva® Channel Separator.

**MATERIAL:** Medical grade plastic.**MARKING:** Label on package: UDI, Lot number and Manufacture date.**INTENDED PURPOSE**

The Motiva® Channel Separator is a surgical instrument, for transient use, intended to create a surgical channel by dissecting the soft tissue from the axilla to the opposing border of the desired breast pocket location to subsequently introduce a sterile breast implant during a transaxillary breast implant procedure.

**INTENDED USERS**

The Motiva® Channel Separator is used in plastic surgeries by qualified surgeons. Not used by lay persons. Not used by patients. Remove the Sterile barrier in the procedure room or operation Theatre.

| CONTRAINDICATION                                 |                                                                                                                |
|--------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| Use for intended purpose only.                   | Use the device for the intended use only. Incorrect use could cause breakage of this device or patient injury. |
| Body area of use.                                | Do not use near central circulatory system and central nervous system.                                         |
| Do not re-use nor re-sterilize.                  | Intended for single use, after use, properly discard according to clinical/practical practice.                 |
| Do not use if package is damaged or compromised. | Sterility could be affected.                                                                                   |
| Do not use after expiration date.                | Sterility may be affected if used after the expiration date stated in the label.                               |

**CAUTIONS/ WARNINGS**

- Before using this product, inspect the package is not damaged or compromised since sterility is affected.
- Keep these instructions accessible to all staff.
- The device should only be used by medical professionals with adequate training in plastic surgery.
- Do not use the system in any other setting than a sterile or environment where proper sterile technique can be applied.
- Do not introduce the handle of the device into the incision site.
- If needed, move the breast tissue upward to shorten the distance between the incision and the end of the channel, avoiding the handle of the device to be pressured against the incision.

**PRECAUTIONS**

- DO NOT modify the Motiva® Channel Separator in any way.
- DO NOT re-use nor re-sterilize, due to the risk of not achieving the expected sterility assurance level.

- Preliminary product examination: Make sure the device and its package is not damaged before using. If these or other defects are present, the device must NOT be used. The device must then be disposed in accordance with clinical/ practice.

**STORAGE**

Store at room temperature (5 to 40°C). Avoid locations with dust, insects and among other vermin. Note that high humidity areas may cause compromising the sterile package in long term storage.

**TRANSPORTATION**

To prevent damage to the medical device during transit, it is recommended the use of racks, trays or rigid containers that protects the device.

**Service life of this product:** Single use, see expiration date on the label.

**CONNECTION TO OTHER DEVICES AND ACCESSORIES**

No other devices are needed to be connected nor related to the Motiva® Channel Separator, the device can be used totally independent.

**RISK WITH THE USE OF THE DEVICE**

No residual risks, nor undesirable side-effects, have been noticed by the manufacturer

|                                          |                                                                                                                                                                                                                                                                                                               |
|------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Physical characteristics</b>          | - The product consists of a single use separator.<br>- Unique mechanical strength resistant plastic material<br>- Specially designed tip and shaft for optimal soft tissue dissection.<br>- Ergonomic and lightweight for improvement surgeon maneuvering performance.<br>- Dimensions: 31cm long X 20mm Wide |
| <b>Operating and Storage Temperature</b> | Temperature: 5 to 40°C                                                                                                                                                                                                                                                                                        |
| <b>Humidity</b>                          | (10 ± 85)% RH                                                                                                                                                                                                                                                                                                 |
| <b>Shelf Life</b>                        | See expiration date on the label.                                                                                                                                                                                                                                                                             |
| <b>Performance Characteristics</b>       | Single use<br>Secure tip                                                                                                                                                                                                                                                                                      |

**GUIDELINES TO CHECK PROPER FUNCTIONALITY OF THE MOTIVA® CHANNEL SEPARATOR**

Check prior to use. Confirm the package is undamaged and its expiration date. Visual check of among others: Straight and not bent, no sharp edges, no splints, no discoloration. Discard damaged products.

**INSTRUCTIONS FOR USE**

- Using sterile technique, open the sterile pack of the Motiva® Channel Separator in the procedure room or operation theatre.
- Place the tip of the Motiva® Channel Separator at the transaxillary incision and enter the incision following the marked path.
- Push the tool until the shaft of the device is located inside the pre-defined opposing border of the marking, and the device's tip reaches the Final Dissection Mark, about 2 cm beneath the edge of the marking. Refers to figure 2.
- Pull out the Motiva® Channel Separator.
- Use precautions according to hospital protocol for the operation and disposal of contaminated waste.


Figure 2. Shaft position

**CLINICAL BENEFITS**

The device does not possess direct clinical benefits. Instead, the device possesses an indirect clinical benefit; assisting plastic surgeon in breast pocket dissection.

**PATIENT POPULATION**

Is intended to be used in women of at least 18 years old undergoing breast implant transaxillary procedure.

**PACKAGING**

1 piece per pack

| LABELLING SYMBOL |                              |                   |                                                                                                 |
|------------------|------------------------------|-------------------|-------------------------------------------------------------------------------------------------|
| <b>UDI</b>       | Unique Device Identification | <b>STERILE EO</b> | Single sterile barrier system with protective packaging outside sterilized using ethylene oxide |
| <b>MD</b>        | Medical Device               |                   | Do not resterilize                                                                              |
| <b>REF</b>       | Catalog Number               |                   | Do not reuse                                                                                    |
|                  | Consult Instructions for use |                   | Temperature Limit                                                                               |
| <b>LOT</b>       | Batch Code                   |                   | Humidity Limitation                                                                             |
|                  | Manufacturer                 |                   | Keep Dry                                                                                        |
|                  | Use by date                  |                   | Keep away from sunlight                                                                         |
|                  | Caution                      |                   | Do not use if package is damaged                                                                |

**REPORTING OF INCIDENT TO MANUFACTURER (SEE BELOW) & COMPETENT AUTHORITY**

In case any patient/user faces a serious incident, please IFU report the incident to the manufacturer, and the Competent Authority of the country where the user/patient resides.

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Copyright© 2024 Establishment Labs S.A.
PulseLavage AB