	IFU Motiva® Reusable Channel Dissector	Document TF04.22.01PLIFUMRCD
	Date Dec 2024	Rev I

EN

Motiva® Reusable Channel Dissector Reprocessing Instructions

Basic UDI	GTIN	PACKING VARIANT	Ref
735009365MRCDDM	07350093650655	Motiva ® Reusable Channel Dissector 20 mm	MRCD-20
735009365MRCDDM	07350093650662	Motiva ® Reusable Channel Dissector 23 mm	MRCD-23


Initial treatment at the point of use: Before using this product, inspect, wash, and sterilize in accordance with these instructions.

Material: Medical grade stainless steel.

Marking: Direct Mark Laser UDI, serial number and Manufacture date.

Intended Use

The Motiva® Reusable Channel Dissector is a reusable surgical instrument, for transient use, intended to create a surgical channel by dissecting the soft tissue from the **transaxillary or inframammary** incision to the opposing border of the desired breast pocket location to subsequently introduce a sterile breast implant during a breast implant procedure.

Contraindication

Use for intended purpose only.	Use the device for the intended use only. Incorrect use could cause breakage of this device or patient injury.
Body area of use.	Do not use near central circulatory system and central nervous system.
Do not use chemicals.	Do not apply chemicals. This could cause corrosion.
Handle with care	Harsh handle can cause deformities or damage the device and will consequently reduce service life.
Do not use polishing powder or wire wool	Do not polish the product's surfaces as this could cause uneven surfaces with sharp edges or it could result in rust or corrosion.
Do not use acid or household detergents	Use only medical neutral detergents (pH 6 to 8) to clean this product. Improper detergent could result in discoloration or corrosion.
Do not use low temperature hydrogen peroxide gas plasma sterilization.	This device is not suitable for sterilization with hydrogen peroxide low temperature gas plasma. It may discolor the surface of the device or affect the feature of the product.

Cautions / Warnings

- Before using this product, inspect, wash, and sterilize in accordance with instructions.
- The success of the reprocessing must be verified by the user.
- Follow instructions and warnings on the cleaning equipment used for reprocessing of this product.
- The quality of the water used for reprocessing shall be compliant to Directive 98/83/EC.

	IFU Motiva® Reusable Channel Dissector	Document TF04.22.01PLIFUMRCD
	Date Dec 2024	Rev I

EN

5. Do not exceed 140°C through reprocessing
6. Keep these instructions accessible to all staff.
7. Do not introduce the handle of the device into the incision site.
8. if needed, move the breast tissue upward to shorten the distance between the incision and the end of the channel, avoiding the handle of the device to be pressured against the incision.


Figure: Shaft position. The tip of the device should be located about 2 cm beneath the edge of the marking.

PRECAUTIONS

1. DO NOT modify the Motiva ® Reusable Channel Dissector in any way.
2. DO NOT use abrasive cleaners or solvents on the device.
3. Preliminary product examination: Make sure the device is not damaged before using. If any defects are present, the instruments must NOT be used. The product must then be disposed in accordance with clinical/practical practice.

Storage

Store in room temperature (5 to 40°C). Avoid locations with dust, insects and among other vermin. Note that high humidity areas may cause corrosion in long term storage.

Transportation

To prevent damage to the medical device during transit, it is recommended the use of racks, trays or rigid containers that protects the device.

Service life of this product: 3 years or 30 reprocessing cycles.

(Subject to following manufacturer's specified maintenance, inspection and proper storage requirements.)

Maintenance / Inspection

Check prior to each use

Visual check of among others: Straight and not bent, no sharp edges, no splints, no corrosion, no discoloration, no scratches, no cracks, and wears.

Discard damaged products.

Preparation for Reprocessing

The instruments must not be disassembled or opened for reprocessing.

Cleaning

The following are the recommended instructions for both processes, manual and automated cleaning. The values/parameters of the processes described here were used in the validation performed and are to be adopted accordingly depending on the existing conditions on site. Do not use fixing agents or hot water (>40°C), as this leads to the fixing of residues and can influence the cleaning success.

Manual Pre-Cleaning

1. Pressure-flush the shaft with a jet pistol 10 seconds (static pressure 3.8 bar).
2. Place the instruments in cold water (10-25°C) for 5 min.
3. Brush the instruments with a soft brush until all visible residues are removed.
4. Pressure-flush all openings and hard-to-reach surfaces and bores with a jet pistol gun for 5 seconds (static pressure 3.8 bar).
5. Pressure-flush the cavities with a jet pistol for 10 seconds.
6. Place the instrument in an ultrasonic bath for 10 min. and sonicate the cleaning solution (0.5% neodisher MediClean forte, 40°C).
7. Pressure-flush the cavities with a jet pistol for 10 seconds (static pressure 3.8 bar).

The next step can be done manually or automatically.

Manual Cleaning

1. Pressure-flush the shaft with a jet pistol for 15 seconds (static pressure 3.8 bar).
2. Place the instruments in the cleaning solution (0.5% neodisher MediClean forte, 40°C) for 5 minutes.
3. Brush the instruments with a soft brush until all visible residues have been removed.
4. Pressure-flush all openings and hard-to-reach surfaces and bores with a jet pistol for 5 seconds each (static pressure 3.8 bar).
5. Pressure-flush the shaft with a jet pistol for 15 seconds.
6. Place the instruments in an ultrasonic bath for 9 min. and sonicate them in the cleaning solution (0.5% neodisher MediClean forte 40°C).
7. Pressure-flush the shaft with a jet pistol for 60 seconds (static pressure 3.6 bar). Then rinse the outer surfaces with cold water for 60 seconds to ensure that the cleaning chemicals are removed.

	IFU Motiva® Reusable Channel Dissector	Document TF04.22.01PLIFUMRCD
	Date Dec 2024	Rev I

EN

Automated Cleaning

Step	Time (min.)	Process Step	Reagents	Temp. (°C)
1	2	Pre-cleaning	Tap water	10-25
2	-	Draining	-	-
3	10	Cleaning	Tap water Dosing: 0,7% Neodisher, MediClean Forte	55
4	-	Draining	-	-
5	3	Neutralization (optional)	Deionized water	10-25
6	-	Draining	-	-
7	2	Rinsing	Deionized water	10-25
8	-	Draining	-	-

Desinfection

Automated and Manual disinfection have been validated using Bomix® plus, BODE Chemie 1 % (v/v) for 300 seconds, then rinse the outer surfaces with cold water for 60 seconds to ensure that the cleaning chemicals are removed. Consider the national regulations and the A0 value = 3000.

Drying

Drying of the outside of the instrument by the drying cycle of the washer/disinfector.

If necessary, additional manual drying can be achieved using a lint-free cloth. Dry hollow spaces of instruments with sterile compressed air.

Functional test, maintenance

Visual inspection for cleanliness according to the operating instructions.

If necessary, repeat the reprocessing process until the instrument is visually clean.

Packaging

The instrument is double wrapped in sterilization pouches (Wipak STERIKING flat rolls type R43 / type R44) according to EN ISO 11607-1, last version.

Sterilization

Sterilizer:

Before preparing the device for sterilization, the surface should be carefully lubricated. Grease-free and temperature-resistant silicone products are recommended. This helps to protect the entire instrument surface from mineral deposits which can later lead to functional impairments. Products used for lubrication must be biocompatible and approved for medical devices. Please note that this instrument lubrication should be performed routinely after each cleaning (ultrasound, etc.) and before each sterilization.

Sterilization of the products with fractionated pre-vacuum process considering the respective national requirements.

- 3 pre-vacuum phases with at least 60 millibar pressure.
- Heating up to a sterilization temperature at 134°C for 5 minutes or at 121°C for 20 minutes.
- Shortest holding time: 3min.
- Drying time: at least 20min.

Note: If contamination with prions (CJD) is suspected, reprocess, and dispose of the instrument.

Additional Instructions

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. If the chemicals and machines described here are not available, and if the preparation process as described is not to be carried out, it is the responsibility of the user to verify and/or validate the process to be used accordingly and to perform routine monitoring of the process.

It is the user's responsibility to ensure that the reprocessing procedure process, including resources, materials, and training of personnel, is appropriate to achieve the required results. The state-of-the-art of the device and national laws require compliance with validated processes.

Limitation of Reprocessing and waste disposal

Repeated processing cycles that include ultrasonic, mechanical washing and sterilization have minimal effects on device's life span and function. The life span is mainly affected by wear or damage during surgical use as well as proper care and handling techniques. After the service life of the instrument is over, it should be submitted for professional disposal or recycling. It shall also be necessary to observe all applicable national laws, regulations or disposable guidelines.

Disclaimer

The end-user must record the reprocessing cycles by Serial Number, ensuring that the device is not used more than the maximum allowed number of 30 reprocessing cycles or 3 years.

Point of Use

Remove the excess soils with a disposable cloth or paper wipe promptly after any use. Do not allow saline, blood, body fluids, tissue, or other organic debris to dry on instruments prior to cleaning.

Maintenance and check by agents

For safe use of this instrument, inspect device prior to and after each use. Alternative or no review and/or inspection may cause injury to the patient and/or healthcare worker and may decrease the performance and function of this device. Additionally, it is recommended to schedule a periodic inspection through your local distributor or the manufacturer.

Patient Population

General population and those prescribed for Plastic Breast Augmentation Surgeries.

	IFU Motiva® Reusable Channel Dissector	Document TF04.22.01PLIFUMRCD
	Date Dec 2024	Rev I

EN

Clinical Benefit

The device is a dissector that does not possess direct clinical benefit. Instead, the device possesses indirect clinical benefits; increase patient satisfaction with better breast pocket dissection.

Packing

1 piece per pack

Return items to your distributor

Defective products must have undergone the entire reprocessing process before being returned.

Used Symbols for IFU and Label

Distributor

Date of manufacture

Unique Device Identification

Consult instructions for use

Medical Device

Non Sterile

Serial Number

Catalog Number

Manufacturer

Caution

Contact your local distributor for more information and report any serious incident involving the device to the manufacturer and the competent authority of the Member State in which the event happened.

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PulseLavage AB