Motiva® Channel Separator
Instructions for Use

**MOTIVA® CHANNEL SEPARATOR
INSTRUCTIONS**

Ref:	MCS-20
Description:	Shaft 20 mm wide, Length 310mm
UDI-DI	07350093650624

Figure 1. Motiva® Channel Separator.

MATERIAL: Medical grade plastic.**MARKING:** Label on package: UDI, Lot number and Manufacture date.**INTENDED USE**

The Motiva® Channel Separator is a single-use surgical instrument, for transient use, intended to create a surgical channel by dissecting the soft tissue from the transaxillary or inframammary incision to the opposing border of the desired breast pocket location to subsequently introduce a sterile breast implant during a breast implant procedure.

INTENDED USERS

The Motiva® Channel Separator is used in plastic surgeries by qualified surgeons. Not used by lay persons. Not used by patients. Remove the sterile barrier in the procedure room or Operation Theatre.

CONTRAINDICATION	
Use for intended purpose only.	Use the device for the intended use only. Incorrect use could cause breakage of this device or patient injury.
Body area of use.	Do not use near central circulatory system and central nervous system.
Do not re-use nor re-sterilize.	Intended for single use, after use, properly discard according to clinical/practical practice.
Do not use if package is damaged or compromised.	Sterility could be affected.
Do not use after expiration date.	Sterility may be affected if used after the expiration date stated in the label.

CAUTIONS/ WARNINGS

- Before using this product, inspect the package and verify it is not damaged or compromised as in this case the sterility will be affected.
- Keep these instructions accessible to all staff.
- The device should only be used by medical professionals with adequate training in plastic surgery.
- Do not open/use the system in any other setting than a sterile OR environment where proper sterile technique can be applied.
- Keep the device handle out of the incision and the tip within the external surgical mark. The shaft of the device should be located about 2cm beneath the edge of the marking.
- If needed, move the breast tissue upward to shorten the distance between the incision and the end of the channel, avoiding the handle of the device to be pressured against the incision.

PRECAUTIONS

- DO NOT modify the Motiva® Channel Separator in any way.
- DO NOT re-use nor re-sterilize, due to the risk of not achieving the expected sterility assurance level.

- Preliminary product examination: Make sure the device and its package is not damaged before using. If these or other defects are present, the device must NOT be used. The device must then be disposed in accordance with clinical practice.

STORAGE

Store at room temperature (5 to 40 °C). Avoid locations with dust, insects and among other vermin. Note that high humidity areas may compromise the sterile package in long term storage.

TRANSPORTATION

To prevent damage to the medical device during transit, it is recommended the use of racks, trays or rigid containers that protect the device.

Service life of this product: Single use, see expiration date on the label.

CONNECTION TO OTHER DEVICES AND ACCESSORIES

No other devices are needed to be connected nor related to the Motiva® Channel Separator, the device can be used totally independent.

RISK WITH THE USE OF THE DEVICE

No residual risks, nor undesirable side-effects, have been noticed by the manufacturer

Physical characteristics	- The product consists of a single use separator. - Unique high mechanical strength plastic material. - Specially designed tip and shaft for optimal soft tissue dissection. - Ergonomic and lightweight for improvement surgeon maneuvering performance. - Dimensions: 31cm long X 20mm Wide
Operating and Storage Temperature	Temperature: 5 to 40 °C
Humidity	(10 - 85)% RH
Shelf Life	See expiration date on the label.
Performance Characteristics	Single use Secure tip

GUIDELINES TO CHECK PROPER FUNCTIONALITY OF THE MOTIVA® CHANNEL SEPARATOR

Check prior to use. Confirm the package is undamaged and its expiration date. Visual check of among others: Straight and not bent, no sharp edges, no splints, no discoloration. Discard damaged products.

INSTRUCTIONS FOR USE

- Before beginning, ensure that the external surgical marks have been created to guide the surgeon during the procedure.
- Using sterile technique, open the sterile pack of the Motiva® Channel Separator in the procedure room or operation theatre.
- Place the tip of the Motiva® Channel Separator at the incision and gently move it along the marked path.
- To accurately dissect the breast pocket, make certain that the device is not inserted beyond the previously created external mark. Advance the instrument until the shaft is within the specified opposing border, ensuring that the tip reaches approximately 2cm outside the edge of the surgical marking, as per the surgical plan. Refer to Figure 2 for guidance.
- Pull out the Motiva® Channel Separator.
- Follow hospital protocol for handling and disposing of contaminated waste.


Figure 2. Shaft position. The tip of the device should be located about 2cm beneath the edge of the marking

CLINICAL BENEFITS

The device does not possess direct clinical benefits. Instead, the device possesses an indirect clinical benefit; assisting plastic surgeon in breast pocket dissection.

PATIENT POPULATION

Is intended to be used in women of at least 18 years old undergoing breast implant procedure.

PACKAGING

1 piece per pack

LABELLING SYMBOL			
	Unique Device Identification		Single sterile barrier system with protective packaging outside sterilized using ethylene oxide
	Medical Device		Do not resterilize
	Catalog Number		Do not reuse
	Batch Code		Temperature Limit
	Manufacturer		Humidity Limitation
	Keep Dry		Use by date
	Caution		Keep away from sunlight
	Do not use if package is damaged		

REPORTING OF INCIDENT TO MANUFACTURER (SEE BELOW) & COMPETENT AUTHORITY

In case any patient/user faces a serious incident, please report the incident to the manufacturer, and the Competent Authority of the country where the user/patient resides.

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Copyright© 2025 Establishment Labs S.A.
PulseLavage AB